

Service Report Preventive Maintenance

Instrument ID: Inova 300
S/N:
Engineer: Don Frank

Date: 3/2/2020
WO#:
OTI No.: 200302001
Probe file: asw cii replace

System Checks:

☐ Clean Console
☒ Pneumatic Operation
☒ Cooling Fans
☐ VT Operation
☐ Clean Filters
ok LHe Boil Off
ok LN Boil Off

Accessory Checks:

☒ PFG Operation - Profile, Recovery, Gshim
☒ Adjust Sample Present on MSR Board
☒ Clean Upper Barrel
☐ SMS - Clean and inspect
☐ Other Accessory - Clean and Inspect

Computer Checks:

CentOS 6	UNIX Version
VnmrJ 4.2A	VNMR Version
ok	Free Disk Space

Other Misc Checks:

☒ Update Shims (/vnmr/shims/*probe*_lk_d2o)
☒ Z0 in probe file (/vnmr/probes/*probe*)
☐ Run Calibrate Macro (C13 Decp, C/H ratio)
☒ Gradient Shim Map
☒ Tune Probe using C6D6

NMR Checks:

	H1	F19	C13	P31	C13 Indirect	Other
PW90	12.5us @ 57		18us @ 62			
S/N	113:1		91:1			
Lineshape	3.13/40.2/69.3					
gamaH2			2502 @ 42			
Bird Test	PASS		PASS			
Grad Strength	20.5 gauss/cm					
GCAL	0.0099903					
GZRecovery	< 1 ms					
Waltz			PASS			
gC/H Ratio						

Comments:

Data Files saved under /home/vnmr1/opentech/asw_cii_replace/benchmark_030220/
d2o shims updated to /vnmr/shims/*probe*_lk_d2o
Probe file /vnmr/probes/*probe* values updated
ASW Probe has broken quartz coil former, needs to be repaired.
PW 90's for low band and high band are long. Suspect AMT amplifier or transmitters have low power.
Remote Status Box has no display.
Missing Lock Reject filter.

OpenVnmrj

1% H2O in D2O doped with GdCl3

50% = 3.13 Hz

.55% = 40.2 Hz

.11% = 69.3 Hz

Sample Name
Date collected 2020-03-02

Pulse sequence **PROTON**
Solvent **d2o**

Temperature **24**
Spectrometer **Inova400.elevance.com-inova400**
Study owner **vnmr1**
Operator **vnmr1**

SAMPLE

date Mar 2 2020
solvent d2o
file /home/vnmr1/Desktop/
Carbonet/opentech/asw_
cii_replace/benchmark_03
0220/shmd2o_030220.fid

DEC. & VT

dn C13
dof 0
dm nnn
dmm c
dmf 200

PROCESSING

fn not used
math i

ACQUISITION


sfrq 299.767
tn H1
at 4.000
np 16000
sw 2000.0
fb 1000
pw 16.0
pw 16.0
tpwr 57
d1 1.000
tof -96.4
nt 1
ct 1
alock n
gain 15

DISPLAY

sp 402.7
wp 1999.8
vs 117
sc 8
wc 234
hzmm 7.23
is 500.00
rfl -402.5
rtp 0
th 20
ins 1.000
nm cdc ph

FLAGS

il n
in n
dp y
hs nn



OpenVnmrJ

1% H2O in D2O doped with GdCl3
H1 PW90
pw=12.5us tpwr=57


Sample Name
Date collected **2020-03-02**

Pulse sequence **PROTON**
Solvent **d2o**

Temperature **24**
Spectrometer **In**

Study owner **vnmr1**
ova300ator **vnmr1**

SAMPLE		PROCESSING	
date	Mar 2 2020	fit	not used
solvent	d2o	math	i
file	/home/vnmr1/Desktop p/Carbonet/opentech/asw_ cii_replace/benchmark_03 0220/h1pw90_030220.fid	werr	
		wexp	
		wbs	
		wnt	
ACQUISITION		DISPLAY	
sfrq	299.767	sp	403.5
tn	H1	wp	1999.0
at	1.000	vs	139
np	4000	sc	8
sw	2000.0	wc	234
fb	1000	h2mm	7.23
pw	arrayed	is	500.00
pw	arrayed	ri	-402.5
tpwr	57	rtp	0
d1	1.000	th	20
tof	-96.4	ins	1.000
nt	1	ai cdc ph	
ct	1		
alock	n		
gain	15		
FLAGS		ACQUISITION ARRAYS	
il	n	array	pw
in	n	arraydim	16
dp	y	i	pw
hs	nn	1	3.125
		2	6.25
		3	9.375
		4	12.5
		5	15.625
		6	18.75
		7	21.875
		8	25
		9	28.125
		10	31.25
		11	34.375
		12	37.5
		13	40.625
		14	43.75
		15	46.875
		16	50
DEC. & VT			
dn	C13		
dof	0		
dm	nnn		
dmm	c		
dmf	200		



OpenVnmrj

Sample Name
Date collected 2020-03-02

Pulse sequence profile
Solvent d2o

Temperature 24
Spectrometer Inova400.elevance.com-inova400
Study owner vnmr1
Operator vnmr1

1% H2O in D2O doped with GdCl3


profile - echo z projection of gradient calibration

profile - gradient calibration

gcal = 0.0099903

max gradient = 20.4601015 gauss/cm

SAMPLE		dm	nnn
date	Mar 2 2020	dmm	c
solvent	d2o	dmm	200
file	/home/vnmr1/Desktop/p/Carbonet/opentech/asw_cii_replace/benchmark_030220/gzprofile_030220.fid	homo	n
	d	dpwr	undefined
ACQUISITION		GRADIENTS	
sfrq	299.767	gzlv1	512
tn	H1	PROCESSING	
at	0.003	sb	0.001
np	512	sbs	0.000
sw	100000.0	wtfile	
fb	55000	proc	ft
tpwr	57	fn	not used
pw	25.0	werr	
p1	12.5	wexp	
d1	1.000	wbs	
d2	0.020	wnt	
tof	-96.4	lp	87409.2
nt	1	rp	-32746.1
ct	1	rp1	0
gain	15	lp1	0
spin	not used	DISPLAY	
FLAGS		sp	-48206.9
il	n	wp	99609.4
in	n	vs	28920
dp	y	sc	8
hs	nn	wc	234
DEC. & VT		hzmm	7.23
dn	C13	is	500.00
dof	0	rfl	48597.5
		rfp	0
		th	23
		ins	1.000
		ai	av



OpenVnmrj

Sample Name
Date collected 2020-03-02


Pulse sequence g2pul
Solvent d2o

Temperature 24
Spectrometer Inova400.elevance.com-inova400

Study owner vnmr1
Operator vnmr1

1% H2O in D2O doped with GdCl3

g2pul - gradient recovery test




OpenVnmrj

1H SENSITIVITY
0.1%Ethylbenzene in CDCl3
S/N = 113:1

Sample Name
Date collected 2020-03-02

Pulse sequence s2pul
Solvent cdcl3

Temperature 24
Spectrometer Inova400.elevance.com-inova400
Study owner vnmr1
Operator vnmr1



OpenVnmrj

Sample Name
Date collected 2020-03-02

Pulse sequence s2pul
Solvent cdcl3

Temperature 24
Spectrometer Inova400.elevance.com-inova400
Study owner vnmr1
Operator vnmr1

1H Bird Test

SAMPLE

date Mar 2 2020
solvent cdcl3
file /home/vnmr1/Desktop/Carbonet/opentech/asw_cii_replace/benchmark_03_0220/h1bird_030220.fid

dmm c
dmf 200

PROCESSING

lb 1.00
fn 32768
math i
werr
wexp
wbs
wnt

ACQUISITION

sfrq 299.761
tn 12.000
at 12.000
np 16026
sw 4006.4
fb 2000
pw 0
tpwr 57
d1 0
tof arrayed
nt 1
ct 1
alock n
gain 40

sp -5504.1
wp 4006.2
vs 119649
sc 8
wc 234
hzmm 3.88
is 500.00
rfi 5504.3
rfp 0
th 20
ins 1.000
ai cdc av

FLAGS


il
in
dp y
hs nn

ACQUISITION APPRAYS

array to
arraydim 3
i tof
1 -5000
2 0
3 5000

DEC. & VT

dn H1
dof -776.0
dm nnn




OpenVnmrj

13C PW90
40% Dioxane in Benzene-d6
pw=18us tpwr=62

SAMPLE	
date	Mar 2 2020
solvent	c6d6
file	/home/vnmr1/Desktop/Carbonet/opentech/asw_cii_replace/benchmark_030220/c13pw90_030220.fid
ACQUISITION	
sfrq	75.381
tn	C13
at	0.400
np	12820
sw	16025.6
fb	9000
pw	18
tpwr	62
d1	5.000
tof	-2113.7
nt	1
ct	1
alock	n
gain	25
FLAGS	
il	n
in	n
dp	y
hs	nn
DEC. & VT	
dn	H1
dof	-121.0
dm	nnn
dmm	c
dmf	200
PROCESSING	
lb	3.50

fn	32768
math	i
werr	
wexp	
wbs	
wnt	
DISPLAY	
sp	2488.7
vp	4671.5
ts	72407
sc	8
wc	234
azmm	18.59
ts	100.00
ff	2983.0
fp	0
ai cdd ph	
ACQUISITION ARRAYS	
array	pw
arraydim	16
i	pw
1	4.5
2	9
3	13.5
4	18
5	22.5
6	27
7	31.5
8	36
9	40.5
10	45
11	49.5
12	54
13	58.5
14	63
15	67.5
16	72

Sample Name	Pulse sequence	Temperature	Study owner
Date collected	Solvent	Spectrometer	Processor
2020-03-02	c6d6	24	vnmr1
		Inova400.elevance.com-inova400	vnmr1



OpenVnmrj

13C SENSITIVITY
40% Dioxane in Benzene-d6
S/N = 91:1

Sample Name
Date collected 2020-03-02

Pulse sequence s2pul
Solvent c6d6

Temperature 24
Spectrometer Inova400.elevance.com-inova400
Study owner vnmr1
Operator vnmr1

SAMPLE

date Mar 2 2020
solvent c6d6
file /home/vnmr1/Desktop/
p/Carbonet/opentech/asw_
cii_replace/benchmark_03
0220/c13sn_030220.fid

ACQUISITION

sfrq 75.383
tn C13
at 0.400
np 12820
sw 16025.6
fb 9000
pw 18.0
pw 18.0
tpwr 62
d1 5.000
tof 0
nt 1
ct 1
alock n
gain 40

FLAGS

il n
in n
dp y
hs nn

DEC. & VT


dn H1
dof -121.0
dm nnn
dmm c
dmf 200

PROCESSING

lb 3.50
fn 32768
math i

DISPLAY

sp -842.9
wp 15999.2
vs 117
sc 8
wc 234
hzmm 40.00
is 100.00
rfl 869.3
rtp 0
th 20
ins 1.000
nm cdc ph



OpenVnmrj

Sample Name
Date collected 2020-03-02

Pulse sequence s2pul
Solvent c6d6

Temperature 24
Spectrometer Inova400.elevance.com-inova400
Study owner vnmr1
Operator vnmr1

¹³C GamaH2


amma-H2 calculation

See J. Magnetic Resonance 7:442 (1972)

J0 = 142.00 Hz, dof[1] = -2003.00 Hz, dof[2] = 997.00 Hz
Jr[1] = 73.00 Hz, Jr[2] = 73.00 Hz

Gamma-H2 = 2502.72 Hz
90 degree pulse = 99.9 usec
coalescence frequency = -503.00 Hz

SAMPLE		dm	y
		dmm	c
date Mar 2 2020		dmf	200
solvent c6d6		dpwr	42
file /home/vnmr1/Desktop/Carbonet/opentech/asw_cii_replace/benchmark_030220/gamah2_030220.fid		PROCESSING	
		lb	3.50
		fn	32768
		math	i
		werr	
		wexp	
		wbs	
		wnt	
ACQUISITION		DISPLAY	
sfrq	75.383	sp	4641.5
tn	C13	wp	829.5
at	0.400	vs	197
np	12820	sc	8
sw	16025.6	wc	234
fb	9000	hzmm	3.54
pw	18.0	is	100.00
pw	18.0	rfl	869.3
tpwr	62	rtp	0
d1	5.000	th	20
tof	0	ins	1.000
nt	1	nm cdc ph	
ct	1	ACQUISITION ARRAYS	
alock	n	array	dof
gain	50	arraydim	2
FLAGS		i	dof
il	n	1	-2003
in	n	2	997
dp	y		
hs	nn		
DEC. & VT			
dn	At		
dof	arrayed		




OpenVnmrj

Sample Name
Date collected 2020-03-02

Pulse sequence s2pul
Solvent c6d6

Temperature 24
Spectrometer Inova400.elevance.com-inova400
Study owner vnmr1
Operator vnmr1

13C Waltz Test



OpenVnmrj

Sample Name
Date collected 2020-03-02

Pulse sequence s2pul
Solvent c6d6

Temperature 24
Spectrometer Inova400.elevance.com-inova400
Study owner vnmr1
Operator vnmr1

13C Bird Test

