5 Replies Latest reply: Jan 14, 2014 2:12 PM by Khan Hekmatyar





Khan Hekmatyar Jan 20, 2014 7:43 PM

Black screen with error message "svc-run returned" - how to fix?

This question has been Answered.

Hello All

It has been running fine till yesterday, I tried to switch to another user account on the vnmrj4 (RHEL 6.3). It rebooted and saying "scv-run returned" in the black screen. it hungs up with the message. I rebooted several times, it is not going anywhere.

Any suggestions.

Thanks

Khan



Correct Answer

by Khan Hekmatyar on Jan 14, 2014 2:12 PM

To fix this,

ssh to your machine after it hangs.

mv /etc/X11/xorg.conf /etc/X11/old_xorg.conf

then rebooted and fixed this issue. It happened due to Nvidia card issue than vnmrJ



Helpful Answer by Khan Hekmatyar

52 Views Categories: Administration: Computer/Printer/Network

Tags: rhel6.3, nvidia, rhel6.x, persistence, video_driver, computer_crash, black_screen, svc-run_returned, video_card, nvidia_card



Josh Kurutz Jan 14, 2014 10:51 AM (in response to Khan Hekmatyar)

Re: svc run

Khan, his sounds distressing. What kind of computer is this? Do you have a CD from which you can boot?





Khan Hekmatyar Jan 14, 2014 11:17 AM (in response to Josh Kurutz)

Re: svc_run

Yes, I have been spending lot of time in fixing one problem to another. it is acquisition computer running red hat 6.3,. No, I do not have CD to boot.

Thanks

Khan





Khan Hekmatyar Jan 14, 2014 2:02 PM (in response to Khan Hekmatyar)

Helpful Answer Re: svc_run

Some how, xorg.conf was corrupted, creating new one will solved this issue

Issue is closed.





Josh Kurutz Jan 14, 2014 2:07 PM (in response to Khan Hekmatyar)

Re: svc_run

Good to hear it's fixed!

How did you create a new xorg.conf file?





Khan Hekmatyar Jan 14, 2014 2:12 PM (in response to Khan Hekmatyar)

Correct Answer Re: svc_run

To fix this,

ssh to your machine after it hangs.

mv /etc/X11/xorg.conf /etc/X11/old_xorg.conf

then rebooted and fixed this issue.It happened due to Nvidia card issue than vnmrJ
Like (0)

Home | Top of page | About Jive | Help

© Agilent 2000-2013 | Powered by jiV2